

Data Sheet -

Pro-Shield Safety Insulin Pen Needle



Product Codes:

AMTL-3104S AMTL-3008S

AMTL-3005S

Latex Free ✓

Phthalates Free ✓

DEHP Free ✓

Device Specifications

Device dimensions (includes cover)	Height	44mm
	Width	16mm
	Length	16mm
Weight (includes cover)		3.17g
Material	Body/Cap/Front shield/Needle hub	Polypropylene
	Lock indicator	ABS
	Sealing film	Medical dialysis paper
	Needle/Springs	Stainless steel
Pen needle	AMTL-3104S	31G 4mm needle length
	AMTL-3005S	30G 5mm needle length
	AMTL-3008S	30G 8mm needle length
Sterilisation method		Ethylene oxide
Intended use		Insulin injection
Connection		Standard universal screw thread
Injection site		(see user instructions)
Product carbon footprint	Functional unit: <i>per pack</i>	TBC gC02e

Packaging Specifications

Pcs/pack		100
Pack weight		356g
Pack dimensions	Height	130mm
	Width	99mm
	Length	135mm
Material		Cardboard (please recycle)
Packs per case		20
Pcs/case		2000
Case weight		7.7kg
Case dimensions	Height	280mm
	Width	280mm
	Length	515mm
Material		Cardboard (please recycle)