

Apollo Pen Needle

Range Compatibility

Apollo BV, Pro-Shield & Pro-Shield Duo pen needles are a universal screw fit.

This compatibility list is not exhaustive.



Novo Nordisk®

NovoPen® 6

NovoPen® Echo Plus

FlexTouch®

FlexPen®

NovoPen® 5

NovoPen® 4

NovoPen® 3

NovoPen® Echo

NovoTwist®

Victoza®

Levemir® FlexPen®

Novolog® FlexPen®

Novolog® Mix 70/30®

NovoPen® Junior

Owen Munford Ltd

Autopen® Classic

Autopen® 24

Sanofi

AllStar® Pro

JuniorSTAR®

ClikSTAR®

AllStar®

Apidra® SoloStar®

Lantus® SoloStar®

Ypsomed AG

Ypsopen®

Eli Lilly & Company Ltd

HumaPen® Savvio™

HumaPen® Luxura™

HumaPen® Luxura™ HD

HumaPen® Ergo II

KwikPen™

Humalog® Mix 50/50

Humalog® Mix 75/25

HumaPen® Memoir™

Humulin® N Pen

Humalin® 70/30 Pen